
*...going one step further*

**P60**

(1012792)

# Simulator for ultrasound guided sclerotherapy for varicose veins

Ultrasound guided foam sclerotherapy for varicose veins is gaining in importance. In order to achieve the optimum result, it is essential that the person carrying out the procedure be suitably qualified. Professional expertise and experience, as well as good hand-eye coordination, are needed to be able to insert the needle accurately into the desired vein. The simulator is a reproduction of the right upper thigh. Its features are:

- It has the natural feel of real soft tissue, ideal for carrying out palpations and injections.
- It is possible to inject real medication
- It features a texture and echogenicity similar to those of real tissue, on realistic ultrasound images
- The simulator is developed and manufactured in Germany

## 1 The 1012792 simulator comprises

A skin-colored right thigh

The following are also delivered

- White tray and cover
- White liquid receptacle
- Stand with plastic bottle and infusion set
- Syringe (5 ml) with cannula (green, 21Gx2)

## 2 Assembly

In order to have access to a vein carrying blood, the liquid supply is to be set up as follows:

- Place the stand on a flat surface
- Fix the knurled screw to hang up the liquid bottle on the upper edge of the stand (marked in red)
- Put 250 ml of water into the plastic bottle and add the red dye until the desired color is achieved for the artificial blood
- Close the bottle with the gray rubber bung, and shake the liquid to produce an even color
- Pull the suspension net over the infusion bottle, and hang it on the hook on the stand
- Pick up the infusion set (4000161) and close the flow regulator
- Insert the spike of the drip chamber right into the gray rubber bung
- Fill it up with liquid from the bottle by pressing on the drip chamber several times, then open the air valve in front of the spike
- Remove the stopper on the end of the infusion cable and attach it to the distal end of the vessel
- Guide the proximal end of the vessel magna into the white receptacle
- Remove any air from the tubes until the artificial blood (red liquid in the plastic bottle) enters the receptacle by opening the flow regulator in the infusion tube
- Close the infusion system once more with the flow regulator
- During the simulation, the flow regulator may be opened for a short time, if necessary, so that the liquid removed may flow
- You can now start the simulation of ultrasound guided sclerotherapy (see chapter 3.4)

## 3 Carrying out ultrasound guided sclerotherapy

### 3.1 Instructions

Vein atrophy (sclerotherapy) is a treatment method whereby varicose veins are „turned off“ by injecting a sclerosing agent, into the vena saphena magna for example. Sclerosing agents lead to a marked modification of the vascular wall. After that, the connective tissue of the veins is modified, and finally becomes atrophied.

Ultrasound guided sclerotherapy is a form of treatment that is being increasingly used, especially in the outpatients' department, by phlebologists. Compared with surgery on varicose veins, sclerotherapy is highly beneficial for the patient in that it is a much less invasive procedure.

The ultrasound guided sclerotherapy method was made public for the first at the conference of the International Union of Phlebology (IUP) in Strasbourg in 1989 [1] and has been described in many published studies as an effective treatment method that has few adverse effects [2,3].

### 3.2 Learning to carry out sclerotherapy on the simulator

In order to learn how to carry out ultrasound guided sclerotherapy and to gain practical experience, dexterity and confidence, it is necessary for the student to get the opportunity to carry out a sufficiently large number of treatments during training.

During initial training, in the interest of patient safety, the necessary techniques and skills can be acquired first of all using a simulator. The simulator must therefore realistically reproduce the circumstances of clinical practice, including anatomy, and the feel and echogenicity of the human skin tissue.

Another advantage of the simulator is that it offers the opportunity to continually repeat the procedure, which in turn increases patient safety. A simulator to learn how to carry out ultrasound guided sclerotherapy is seen by experienced practicing phlebologists as a useful pedagogical tool [4].

### 3.3 What is needed for simulation?

To simulate ultrasound guided sclerotherapy, you need, in addition to the prepared simulator shown here:

- An ultrasound device with a 10-18 MHz probe
- Materials to prepare and carry out sclerotherapy (syringe, cannula and medication or a suitable substitute)

The innovative material that the simulator is made of is unusually realistic thanks to:

- The consistency of the material
- The material's resistance
- The echogenicity of the material

#### Consistency of the material

The consistency of the material while positioning the probe and injecting the material guided by ultrasound, corresponds to the consistency of human soft tissue, so for the student, it feels like carrying out a real procedure on a patient.

#### The material's resistance

The resistance of the material used for ultrasound guided puncture of the vessel replicates the features of real human tissue perfectly.

#### Echogenicity of the material

The ultrasound characteristics of the material on the simulator are very similar to those of a real life situation on a patient.

### 3.4 Carrying out the simulation step by step

Carrying out the simulation of ultrasound guided sclerotherapy

- Connect the accessories with the simulator as described below: The "blood" should flow – just like in the veins – in a distal to proximal direction, so the following position should be adopted by the doctor/user:

For right-handed people:

Infusion bottle with the infusion system on the right, and the reservoir on the left (the user is placed in a dorsal position)

For left-handed people:

Infusion bottle with the infusion system on the left, and the reservoir on the right (the user is placed in a ventral position)

Fig. 1

Simulator for ultrasound guided sclerotherapy

## English

- Fill up the infusion bottle with artificial blood and then remove all air from the infusion system and from the artificial vessel system of the simulator until the artificial blood starts to flow into the reservoir (see also chapter 2)

Fig. 2

Simulator after set-up with a full infusion and vessel system

- Sweep over the surface of the simulator with the probe in the area where the cannula is to be inserted.
- Fill the syringe with the desired sclerosing agent and put it down ready for use. Use an injection cannula of between 30 and 50 mm in length.
- Position the ultrasound probe and try to find the longitudinal direction in which the puncture is to be carried out.

Fig. 3 [4]

A longitudinal section of the ultrasound image of the artificial vena saphena magna

- To keep the ultrasound probe stable, support the hand, and hold in the other hand the syringe that has been filled with the sclerosing agent.
- Carry out the puncture on the upper thigh with the cannula, approximately in the middle of the ultrasound probe, at an angle of 45°.

Fig. 4 [4]

Carry out the ultrasound guided injection on the simulator with the syringe at a 45° angle

- Watch the path of the cannula on the screen of the ultrasound device, and carefully move the cannula longitudinally towards the vessel wall.
- As soon as the vessel wall has been punctured, move the cannula forwards until the tip and bevel can be seen in the middle of vascular lumen.

Fig. 5 [4]

Ultrasound image of the correct intraluminal position of the tip of the cannula in the artificial vena saphena magna

- Aspirate the artificial blood visibly into the syringe to make sure that the cannula has been placed in the correct intraluminal position.

Fig. 6 [4]

Aspirate the artificial blood into the syringe before injecting the sclerotherapy medication

- Hold the position of the ultrasound probe and the cannula, and follow its position on the screen of the ultrasound device. Inject the sclerosing agent slowly and continuously. If you can make out a white striping (hyperdense = highly echogenic, depending on the medication used), on the screen of the ultrasound device, in the vessel lumen, you have successfully completed the ultrasound guided sclerotherapy simulation.

Fig. 7 [4]

Ultrasound guided injection of the medication longitudinally along the vena saphena magna on the simulator

### 3.5 Summary

The simulator to learn how to carry out ultrasound guided sclerotherapy is suitable both for the purposes of training students and for advanced training. During training, thanks to the simulator, the student will become familiarized with the basic movements and the procedure, so that inexperienced practitioners will not have to learn by carrying out their first procedure on a patient, thus avoiding the possibility of causing harm.

For advanced training, the simulator provides an opportunity for practicing phlebologists and vascular surgeons to hone their techniques, and to gain confidence, through professional development, in carrying out ultrasound guided sclerotherapy, in order to minimize the risk of intra-arterial or extra-arterial injection on a real patient.

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### 4 Cleaning, taking care of and storing the simulator

After use, the tubes, the thigh and the receptacle are to be cleaned.

They are to be cleaned first of all under running water:

- The plastic bottle
- The thigh
- The receptacle

To clean, at most use a little washing up liquid instead of just pure water.

Then fill the bottle once more with water, and add a drop of washing up liquid. Assemble the infusion set, simulator and receptacle, and let the liquid flow through the entire system once.

Then put the parts in a well-ventilated place, out of direct sunlight, for them to dry. Dry the outside of the thigh, rub in talcum powder (shaking off any excess), and place on the white tray (so that no pressure marks form). Protect the thigh with the cover supplied for storage.

### 5 Product information

Box item delivered (16x60x34 cm)

Size: L46 X H12 X W12.5 cm (with cover)

Weight: 5.4 kg.

Storage temperature: < 50°C

Working temperature: 20-40°C

Do not store the device in direct sunlight.

### 6 Guarantee and Certifications

3 year guarantee not valid for consumables and worn parts.

### 7 Accessories and replacement parts (Fig.1)

Product	Article n°	Quantity in the set
① Artificial blood concentrate, 10 ml	4000160	10 ml
② Plastic bottle with mounting, 250 ml	4000155	1
③ Infusion set	4000161	1
④ Gray rubber bung	1012777	1
⑤ Tube connector	1013590	20
⑥ Tigh (silicone mix) with pipe adapter (2),PS protective cover	1013594	1

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© Copyright 2012 for instruction manual and design of product:
**3B Scientific GmbH**
3B Scientific GmbH, Germany